



## User's Manual

# p62 ELISA Kit



DEIA6457



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The p62 Enzyme-Linked Immunosorbent Assay (ELISA) kit is a complete kit for the quantitative determination of p62 in cell lysates of human, mouse and rat origin.

### Storage

All components of this kit, except the Standard, are stable at 4°C until the kit's expiration date. The standard must be stored at or below -70°C upon receipt. For more detailed information, please download the following document on our website.

### Sensitivity

100 pg/mL