



User's Manual

## SMN ELISA Kit



DEIA6483



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The SMN ELISA kit is a complete kit for the quantitative determination of SMN in samples of human and mouse origin.

### Storage

All components of this kit, except the Standards, are stable at 4°C until the kit's expiration date. The Standards should be stored at or below -20°C upon receipt. For more detailed information, please download the following document on our website.

### Sensitivity

50 pg/mL

