



User's Manual

Horse IL-2 ELISA Kit



DEIA649



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Enzyme linked immunosorbent Assay(ELISA) for the detection of Horse IL-2 in cell culture supernatants. Kit contains sufficient components to quantitate Horse IL-2 protein concentration in up to 40 samples, tested in duplicate.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

15-0.23 ng/mL

