



User's Manual

Human sTNF-R (60kDa) Platinum ELISA Kit



DEIA6494



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human sTNF-R (60kDa) ELISA is an enzyme-linked immunosorbent assay for the quantitative detection of Human sTNF-R (60kDa).

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.