



User's Manual

arg8-Vasopressin ELISA Kit



DEIA6498



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Arg8-Vasopressin EIA kit is a competitive immunoassay for the quantitative determination of Vasopressin in samples.

Storage

All components of this kit, except the conjugate and standard, are stable at 4°C until the kit's expiration date. The conjugate and standard must be stored frozen at -20°C. For more detailed information, please download the following document on our website.

Sensitivity

3.39 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com