



User's Manual

Human APOH ELISA Kit



DEIA685



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This enzyme linked immunosorbent assay (ELISA) is for the detection of Human Beta-2 Glycoprotein 1 (β 2 GP1) in serum, plasma, cerebrospinal fluid, milk, and urine. This kit contains sufficient components to quantitate Human β 2GP1 in up to 40 samples, tested in duplicate.

Storage

Store components of the kit at 2-8°C upon arrival up to the expiration date. For more detailed information, please download the following document on our website.



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com