



User's Manual

Human Acrogranin (GRN) ELISA KIT

REF

DEIA6919



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Human Acrogranin (GRN) ELISA KIT is suitable for the quantitative determination of human Acrogranin/Progranulin(PGRN) concentrations in cell culture supernates, serum and plasma.

Storage

Store at 2-8°C for up to 6 months. For more detailed information, please download the following document on our website.

Precision

Intra-assay Precision: 6-8%

Inter-assay Precision: 10-12%

Detection Range

0.32-200 ng/mL

Sensitivity

1.6 ng/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com