



## User's Manual

# Human Cerebellin-1 ELISA Kit



DEIA6949



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Human Cerebellin-1 ELISA has been shown to accurately quantitate the recombinant and natural human Cerebellin-1.

### Storage

Store at 2-8°C for up to 6 months. For longer storage, unopened Standard, Positive Control and Biotin Solution Concentrated should be stored at -20 or -70°C for up to 8 months.

### Precision

Intra-assay Precision: 6-8%

Inter-assay Precision: 8-12%

### Detection Range

Standard range: 0.16-100 ng/mL

Dynamic range: 0.8-100 ng/mL

### Sensitivity

0.8 ng/mL