



User's Manual

Human Amyloid β 40 High Sensitivity ELISA KIT

REF

DEIA7362



96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Human Amyloid β 40 High Sensitivity ELISA KIT is used for the measurement of Amyloid beta 40 in cerebrospinal fluid, cell culture supernatants, primary neurons and plasma.

Storage

Recommended storage for kit components is 2-8°C.

All components are shipped and stored at 2-8°C. Reconstituted standards and quality controls can be frozen for future use, but repeated freeze thaws should be avoided. Refer to expiration dates on all reagents prior to use. Do not mix reagents from different kits unless they have the same lot numbers.

Precision

Intra-assay (% CV): <10%

Inter-assay (% CV): <10%

Detection Range

16-500 pg/mL

Sensitivity

4.0 pg/mL