



**User's Manual**

# Brain Derived Neurotrophic Factor ELISA Kit



**DEIA7373**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

BDNF kit is designed to measure the amount of BDNF in cell culture supernatants, tissue homogenates and biological fluid (serum, plasma, and serum-free) samples from Human and rat. There are enough reagents included in this kit for two 96-well immuno-assay plates.

### Storage

Maintain the unopened kit at 2-8°C until expiration date indicated on the label. For more detailed information, please download the following document on our website.

### Detection Range

7.8 pg/mL to 500 pg/mL

### Sensitivity

7.8 pg/mL

