



User's Manual

Pigment Epithelium Derived Factor ELISA Kit



DEIA7374



96T






This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Pigment Epithelium Derived Factor ELISA KIT is designed to measure the amount of PEDF in cell culture supernatants, tissue homogenates and biological fluid (vitreous, serum, plasma, and serum-free) samples of human origin.

Storage

Maintain the unopened kit at 2-8°C until expiration date indicated on the label. For more detailed information, please download the following document on our website.

Detection Range

0.9-62.5 ng/mL

Sensitivity

0.9 ng/mL

