



User's Manual

Glucagon ELISA Kit



DEIA7396



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Glucagon ELISA KIT is used for the non-radioactive quantification of intact glucagon level in human, rat, mouse, and porcine serum/plasma samples.

Storage

Prior to use, all components in the kit can be stored up to 2 weeks at 2-8°C. For longer storage (>2 weeks), freeze diluted Wash Buffer, Assay Buffer, Standards and Quality Controls at $\leq -20^{\circ}\text{C}$. For more detailed information, please download the following document on our website.

Sensitivity

0.003 ng/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com