



User's Manual

## Ras Activation ELISA Kit



DEIA7400



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Ras Activation ELISA KIT is designed to detect the presence of activated Ras.

### Storage

Kit components arrive on dry ice, and must immediately be stored at the temperatures specified below. Kit components are stable for 3 months from date of shipment when stored as directed. For more detailed information, please download the following document on our website.