



User's Manual

Human Fibroblast Growth Factor-23 (FGF-23) ELISA Kit



DEIA7424



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human Fibroblast Growth Factor-23 (FGF-23) ELISA Kit is intended for the non-radioactive quantification of Human FGF-23 in serum, plasma, and adipocyte extracts or cell culture media samples.

Storage

Store kit and unopened components at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

3.5 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com