



User's Manual

# Mouse Angiopoietin-Like 3 ELISA KIT



DEIA7461



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Mouse Angiopoietin-Like 3 ELISA KIT is to be used for the quantitative determination of mouse/ rat ANGPTL3 in biological fluids.

### Storage

Reagents must be stored at 2-8°C when not in use. Bring reagents to room temperature before use. Do not expose reagents to temperatures greater than 25°C.

### Precision

Intra-assay (CV%): 3.72%-5.88%

Inter-assay (CV%): 5.07%-7.96%