



User's Manual

Entamoeba histolytica IgG ELISA Kit



DEIA7704



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Entamoeba histolytica IgG ELISA KIT is intended for the qualitative determination of antibodies against Entamoeba histolytica antigens in human serum or plasma (citrate).

Storage

The reagents are stable up to the expiry date stated on the label when stored at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is >95%.