



User's Manual

HORSE HAPTOGLOBIN ELISA Kit



DEIA7710



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Horse Haptoglobin ELISA kits are a highly sensitive two- site enzyme linked immunoassay (ELISA) intended for measuring Haptoglobin in horse biological samples.

Storage

Store component of the kit at 2-8°C. For more detailed information, please download the following document on our website.