



User's Manual

Human HGF ELISA Kit



DEIA7945



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human HGF ELISA KIT is suitable for the detection and quantification of human hepatocyte growth factor (hHGF) concentrations in cell lysates, sera and plasma.

Storage

Note: If used frequently, reagents may be stored at 2-8°C. If used infrequently, reagents should be stored at -20°C.

Sealed, Unopened Assay Kit: 2-8°C, 1 month.