



User's Manual

Human MIF ELISA Kit



DEIA8068



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Human MIF ELISA Kit is intended for detection and quantification of Human Macrophage Migration Inhibitory Factor(hMIF) Concentrations in Cell Lysates, Sera and Plasma.

Storage

1. Sealed, Unopened Assay Kit: 2-8°C/1 month For more detailed information, please download the following document on our website.