



## User's Manual

# Human NCALD ELISA KIT



DEIA8714



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Human NCALD ELISA KIT is for the quantitative detection of human NCALD in serum and plasma.

### Storage

Microplate and Standards can be store at -20°C for six months, but the other materials should be store at 2-8°C for six months.

### Precision

Intra-assay (CV%): 3.9%-6.9%

Inter-assay (CV%): 4.5%-5.0%

### Sensitivity

3 pg/mL