



## User's Manual

# Human SERPINF2 ELISA Kit



DEIA9395



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Human SERPINF2 ELISA Kit is for the quantitative determination of active alpha-2-antiplasmin in human plasma.

### Storage

All kit components must be stored at 4°C. Kit should be used no later than the expiration date.

### Detection Range

The assay measures active alpha-2-antiplasmin in the 0-10 ng/mL range. Extensive dilutions must be performed to assure the unknowns will be in the assay range. It is highly suggested to dilute unknowns by at least 1:1,000 because of the high level of alpha-2-antiplasmin in normal human plasma and serum.