



**User's Manual**

# Human Cytochrome C ELISA Kit



**DEIA9427**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Human Cytochrome C ELISA Kit is an enzyme-linked immunosorbent assay for the quantitative detection of human cytochrome c in cell culture lysates, human whole blood or serum.

### Storage

Store kit reagents at 2-8°C. Immediately after use, remaining reagents should be returned to cold storage (2-8°C). Expiry of the kit and reagents is stated on box front labels. The expiry of the kit components can only be guaranteed if the components are stored properly, and if, in case of repeated use of one component, the reagent is not contaminated by the first handling.

### Precision

Intra-assay coefficient of variation: 6.0%

Table 5.

Inter-assay coefficient of variation: 3.96%

Table 6.

### Sensitivity

0.05 ng/mL