



User's Manual

Human sCD95 (APO1/Fas) ELISA KIT



DEIA9438



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human sCD95 (APO1/Fas) ELISA KIT is to be used for the quantitative determination of soluble CD95 (APO-1, Fas) in human serum, plasma, buffered solutions or cell culture medium. The assay will recognize both natural and recombinant human sCD95.

Storage

Store kit at 2-8°C or -20°C upon arrival up to the expiration date. Avoid multiple freeze/thaw cycles. For more detailed information, please download the following document on our website.

Precision

Detection Range

93.75 to 3000 pg/mL

Sensitivity

47 pg/mL