



User's Manual

Human IGFBP-2 ELISA Kit



DEIA9469



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human IGFBP-2 ELISA Kit is an enzyme-linked immunosorbent assay for the quantitative measurement of Human IGFBP2 in serum, plasma, cell culture supernatants and urine.

Storage

Kit may be stored for up to 6 months at 2-8°C from the date of shipment. Standard (recombinant protein) should be stored at -20°C or -80°C (recommended at -80°C) after reconstitution. For more detailed information, please download the following document on our website.

Sensitivity

20 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com