



User's Manual

Human Transforming Growth Factor Alpha ELISA KIT



DEIA9578



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human Transforming Growth Factor Alpha ELISA KIT is to be used for the quantitative measurement of Human TGF alpha in serum, plasma, cell culture supernatants and urine.

Storage

Stable for up to 6 months from date of shipment at 2-4°C. Store reconstituted standard (recombinant protein) at -80°C. Opened Microplate Wells and reagents are stable for 1 month at 2-4°C. Return unused wells to the pouch containing desiccant pack and reseal along the entire edge. For more detailed information, please download the following document on our website.

Sensitivity

3 pg/mL