



User's Manual

Human THPO ELISA Kit



DEIA9584



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human THPO ELISA Kit is an enzyme-linked immunosorbent assay for the quantitative measurement of human TPO in serum, plasma, cell culture supernatants, and urine.

Storage

Stable for up to 6 months from date of shipment at 2-4°C. Store reconstituted standard (recombinant protein) at -80°C. Opened Microplate Wells and reagents are stable for 1 month at 2-4°C. Return unused wells to the pouch containing desiccant pack and reseal along the entire edge.

Note: The whole kit is stable for 1 year when stored at -20°C. Avoid repeated freeze-thaw cycles.

Sensitivity

25 pg/mL