



User's Manual

Monkey IL-23 ELISA kit



DEIA9606



5 x 96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Monkey IL-23 ELISA kit is used for the quantitative detection of monkey IL23 in culture supernatant.

Storage

Store all reagents at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

5-320 pg/mL

Sensitivity

5 pg/mL