



User's Manual

Mouse Interferon alpha ELISA Kit



DEIA9650



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Mouse Interferon alpha ELISA Kit is suitable for the detection of IFN-alpha (IFN- α) in serum and plasma samples.

Storage

For retention of activity, all reagents should be kept at 2-6°C in the dark. Deionized or distilled water should be used for preparation of all reagents. All dilutions should be made with polypropylene tubes and pipette tips. Pipette tips should be changed between each dilution tube. All measurements for standards and samples should be performed in duplicate. At least two control wells (wells with Sample Buffer only) should be used for each assay; these control values should be subtracted from all readings prior to any calculations or plots of the data. For more detailed information, please download the following document on our website.

Detection Range

12.5-400 pg/mL

