



User's Manual

Mouse Prorenin/Renin Total Antigen ELISA Kit



DEIA9664



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics



Address: 45-1 Ramsey Road, Shirley, NY 11967, USA



Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com



Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Mouse Prorenin/Renin Total Antigen ELISA Kit is intended for the quantitative determination of total prorenin and renin antigen in mouse plasma and serum.

Storage

All kit components must be stored at 4°C. Store unopened plate and any unused microtiter strips in the pouch with desiccant. Reconstituted standards and primary may be stored at -70°C for later use..

Detection Range

The assay measures total mouse prorenin and renin in the 0.1-100ng/mL range.