



User's Manual

Cr(Creatinine) Assay Kit(Colorimetric)



DEIABL-M4



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Cr Assay Kit allows for the in vitro quantitative determination of Cr concentrations in serum, plasma, tissue homogenates and other biological fluids.

Storage

2-8 °C for 6 months

Precision

Intra-Assay: CV<8%

Inter-Assay: CV<10%

Detection Range

25-800µmol/L

Sensitivity

< 10µmol/L