



User's Manual

## Sulfonamides ELISA Kit



DEIABL-QB10



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Sulfonamides ELISA Kit(15-in-1) is developed for quantitative measurement of Sulfonamides in tissue, liver (chicken).

### Storage

Store the reagents at 2-8°C until expiration date.

### Detection Limit

Tissue(High detection limit)- 1ppb

### Sensitivity

1ppb