



User's Manual

Oxacillin ELISA Kit

REF DEIABL-QB30

 96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Oxacillin ELISA Kit is developed for quantitative measurement of Oxacillin in tissue, milk.

Storage

Store the reagents at 2-8°C until expiration date.

Detection Limit

Milk-5ppb; Tissue-2ppb

Sensitivity

0.5ppb