



## User's Manual

# Lincomycin ELISA Kit



DEIABL-QB33



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Lincomycin ELISA Kit is developed for quantitative measurement of Lincomycin in tissue, liver, fish, shrimp, honey.

### Storage

Store the reagents at 2-8°C until expiration date.

### Detection Limit

Tissue, Fish, Shrimp, Royal Jelly- 4ppb

### Sensitivity

0.2ppb