



## User's Manual

# Ethoxyquin ELISA Kit

REF

DEIABL-QB36



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Ethoxyquin ELISA Kit is developed for quantitative measurement of Ethoxyquin in tilapia, shrimp.

### Storage

Store the reagents at 2-8°C until expiration date.

### Detection Limit

5ppb

### Sensitivity

1ppb