



User's Manual

Benzimidazoles ELISA Kit

REF

DEIABL-QB43



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Benzimidazoles ELISA Kit is developed for quantitative measurement of Benzimidazoles in tissue, milk, eggs.

Storage

Store the reagents at 2-8°C until expiration date.

Detection Limit

10ppb

Sensitivity

0.5ppb