



User's Manual

Sarafloxacin ELISA Kit



DEIABL-QB5



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Sarafloxacin ELISA Kit is developed for quantitative measurement of Sarafloxacin in fish, shrimp.

Storage

Store the reagents at 2-8°C until expiration date.

Detection Limit

0.4ppb

Sensitivity

0.2ppb