



User's Manual

Carbofuran ELISA Kit

REF

DEIABL-QB57



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Carbofuran ELISA Kit is developed for quantitative measurement of Carbofuran in tea, vegetable.

Storage

Store the reagents at 2-8°C until expiration date.

Detection Limit

Tea: 5ppb; Vegetable: 10ppb

Sensitivity

0.5ppb