



Rabbit Anti-Peg-filgrastim polyclonal antibody (DEIABL218P)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody binds to PEG filgrastim
Target	PEG Filgrastim
Immunogen	PEG filgrastim
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	Lot specific
Size	100 µg
Buffer	Phosphate buffered saline (PBS), BSA, glycerol
Preservative	See individual product datasheet
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Pegfilgrastim is a PEGylated form of the recombinant human granulocyte colony-stimulating factor (GCSF) analog filgrastim. It serves to stimulate the production of white blood cells (neutrophils). Pegfilgrastim treatment can be used to stimulate bone marrow to produce more neutrophils to fight infection in patients undergoing chemotherapy.
Keywords	Granulocyte colony-stimulating factor pegfilgrastim;Peg-filgrastim;pegfilgrastim-cbqv;pegfilgrastim-jmdb

GENE INFORMATION

Synonyms	Granulocyte colony-stimulating factor pegfilgrastim; Peg-filgrastim; pegfilgrastim-cbqv; pegfilgrastim-jmdb
----------	---