



User's Manual

Bovine IL-6 ELISA Kit



DEIABL385



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Interleukin 6 (IL-6) test kit is a highly sensitive twosite enzyme linked immunoassay (ELISA) for measuring IL-6 in biological fluid of Bovines.

Storage

4°C

Detection Range

0.0625 ng/ml - 4 ng/ml