



User's Manual

## Cat Hemoglobin ELISA Kit



DEIABL387



96T



This product is for research use only and is not intended for diagnostic use.


For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Hemoglobin test kit is a highly sensitive two-site enzyme linked immunoassay (ELISA) for measuring hemoglobin in biological fluids of cats.

### Storage

4°C

### Detection Range

12.5 ng/ml - 400 ng/ml