



User's Manual

Chicken Ovotransferrin ELISA Kit



DEIABL390



96T




This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The total Chicken Ovotransferrin (TX) test kits are a highly sensitive two-site enzyme linked immunoassay (ELISA) for measuring Ovotransferrin (Conalbumin) in biological samples of Chickens.

Storage

4°C

Detection Range

6.25 ng/ml - 400 ng/ml