



**User's Manual**

# **Dog Retinol Binding Protein ELISA Kit**



**DEIABL398**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Retinol Binding Protein (RBP) test kit is a highly sensitive two-site enzyme linked immunoassay (ELISA) for measuring RBP in biological fluid of Dogs.

### Storage

4°C

### Detection Range

3.125 ng/ml - 100 ng/ml