



## User's Manual

# Pig Alpha 1-Acid Glycoprotein ELISA Kit

REF

DEIABL412



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The alpha 1-Acid Glycoprotein (AGP) test kit is a highly sensitive two-site enzyme linked immunoassay (ELISA) for measuring AGP in biological fluid of Pigs.

### Storage

4°C

### Detection Range

6.17 ng/ml - 1500 ng/ml