



User's Manual

Rat Alpha 1-Acid Glycoprotein ELISA Kit

REF

DEIABL413



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Alpha 1-Acid Glycoprotein test kits are a highly sensitive two-site enzyme linked immunoassay (ELISA) for measuring Alpha 1-Acid Glycoprotein in biological samples of rats.

Storage

4°C

Detection Range

12.5 ng/ml - 400 ng/ml