



User's Manual

Rat GST Yb1 (mu) ELISA Kit



DEIABL414



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Glutathione S-Transferase Yb1 (GST Yb1) test kit is a highly sensitive two-site enzyme linked immunoassay (ELISA) for measuring Mu glutathione S-transferase (GST Yb1, μ GST) in biological samples of rats.

Storage

4°C

Detection Range

3.125 ng/ml - 100 ng/ml