



User's Manual

Rat Hemopexin ELISA Kit



DEIABL416



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Hemopexin test kits are a highly sensitive twosite enzyme linked immunoassay (ELISA) for measuring Hemopexin in biological samples of rats.

Storage

4°C

Detection Range

6.25 ng/ml - 400 ng/ml