



User's Manual

Sheep Haptoglobin ELISA Kit



DEIABL419



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Haptoglobin (HPT) test kit is a highly sensitive twosite enzyme linked immunoassay (ELISA) for measuring HPT in biological fluid of Sheep.

Storage

4°C

Detection Range

15.625 ng/ml - 1000 ng/ml

