



# Rabbit Anti-E. coli HCP Polyclonal Antibody (DEIABL481D)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	E. coli HCP
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Purified
Conjugate	Unconjugated
Applications	WB, ELISA Western blot detection and characterization of E. coli host cell proteins. This antibody will detect E. Coli HCP in ELISA kit (Catalog: # DEIABL481).
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Short Term: 2°C to 8°C
Ship	Wet ice

## BACKGROUND

<b>Introduction</b>	Host cell proteins (HCPs) are process-related protein impurities that are produced by the host organism during biotherapeutic manufacturing and production.
---------------------	---

---

<b>Keywords</b>	HCP;Host Cell Protein
-----------------	-----------------------

---