



User's Manual

Human Hyaluronidase1/HYAL1 ELISA Kit

REF

DEIABL640



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Sandwich High Sensitivity ELISA kit for Quantitative Detection of Human Hyaluronidase1/HYAL1.

Reagents And Materials Provided

- 96-well plate precoated with anti-Human HYAL1 antibody: 1
- lyophilized recombinant Human HYAL1 standard: 10ng/tube
- biotinylated anti-Human HYAL1 antibody: 130ul (dilution 1:100)
- Avidin-Biotin-Peroxidase Complex(ABC): 130ul (dilution 1:100)
- sample diluent buffer: 30ml
- antibody diluent buffer: 12ml
- ABC diluent buffer: 12ml
- TMB color developing agent: 10ml
- TMB stop solution: 10ml
- Adhesive cover: 4

Storage

Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles(Shipped with wet ice.)

Detection Range

15.6pg/ml-1000pg/ml

Sensitivity

<10pg/ml



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com