



## User's Manual

# Mouse IL-17RC ELISA Kit



DEIABL643



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Sandwich High Sensitivity ELISA kit for Quantitative Detection of Mouse IL-17RC.

### Reagents And Materials Provided

- 96-well plate precoated with anti-Mouse IL17rc antibody: 1
- lyophilized recombinant Mouse IL17rc standard: 10ng/tube
- biotinylated anti-Mouse IL17rc antibody: 130ul (dilution 1:100)
- Avidin-Biotin-Peroxidase Complex(ABC): 130ul (dilution 1:100)
- sample diluent buffer: 30ml
- antibody diluent buffer: 12ml
- ABC diluent buffer: 12ml
- TMB color developing agent: 10ml
- TMB stop solution: 10ml
- Adhesive cover: 4

### Storage

Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles(Shipped with wet ice.)

### Detection Range

156pg/ml-10,000pg/ml

### Sensitivity

<10pg/ml